# 技术指南 Technical Note

PA\_2022001\_00

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### Usage and maintenance of electrodes

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With the continuous improvement of glass formulation and the continuous improvement of the production process, the performance of pH glass electrodes is getting better and better. And because the pH glass electrode is stable compared to other pH electrodes and is not easily interfered with, the pH glass electrode is currently one of the most widely used electrodes.

However, the usage of pH glass electrodes has a certain lifespan, which depends on the nature of the object to be measured, the temperature of the analyte, and the use and maintenance of the electrode. Then the use and maintenance of electrodes is even more important, this article will introduce how to use and maintain the electrodes correctly and extend their service life.

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#### Usage and maintenance of electrodes

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Rinse the electrode with distilled water before and after the measurement, and then use a swab cloth without hemp to suck away the excess water on the electrode head, do not rub, otherwise it will generate static electricity and interfere with the accurate measurement of pH glass film.

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When store the pH electrodes, keep wet, it is recommended to put in the electrode soaking bottle, and there is a special soaking solution in the bottle, which can protect the measuring electrode and the reference electrode, do not put the electrode in distilled water or deionized water storage, which will lead to accelerated failure of the electrode, if in the storage found that there is a white potassium chloride crystal volume on the electrode, the white crystal substance will not affect the measurement, just use deionized water to rinse the electrode to remove the crystal after sucking dry and then use it.

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When using a rechargeable pH electrode, the reference solution should be added to the dosing hole, but not exceeded, and the filling hole plug should be opened during measurement to ensure that the reference solution is easy to flow out of the salt bridge. However, for the special structure of the tris electrode, it is recommended not to open the liquid filling plug when measuring the liquid to be measured, so as to avoid the flow rate being too fast and affecting the measurement accuracy.

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The most common problem in pH measurements is reference salt bridge obstruction, whose symptoms include slow response, measurement out of normal range, noise generation, different types of electrode references, and different maintenance regimens:

a: Gel-filled (non-rechargeable): Immerse the electrode in a beaker filled with 60°C of warm water for 15 min, remove the gel or salt that dried up on the salt bridge, and then place the electrode into a beaker containing 4

mol of KCl warm solution until it returns to room temperature.

b: Liquid-filled type (rechargeable type), drain the electrolyte in the electrode, clean the cavity with distilled water, re-inject the fresh electrolyte, and then soak the electrode in 60 °C fillable liquid for 15 minutes.

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After a long time, if there is dirt on the glass film, it is recommended to clean it according to the following method:

A: If the glass film attachment is protein, it can be cleaned with soap liquid, gently rub the glass film with a soft cloth, pay attention to slight uniform force when wiping, the bubble film is thin and easy to break, pay attention to protection.

b: If it is an inorganic substance, it can be cleaned with 0.1mol hydrochloric acid and then rinsed with deionized water.

c: If it is oily, it can be cleaned with acetone or methanol, do not soak, or can 1000104603

be washed with soap liquid, and then rinsed with deionized water.

After cleaning the electrodes in any of the above schemes, the electrodes should be soaked in the storage soaking solution for more than 15 minutes before use.

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When the electrode is used in the test, pay attention to the protection of the BNC plug to prevent the test solution or other liquid from splashing on the plug, to ensure that the plug is dry, to avoid bad insulation between the BNC plug shell and the inner core, resulting in the impact of test performance. If the liquid is touched, it can be sucked dry in time, and the moisture can also be dried with a hair dryer.

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##### Proper usage and storage of pH buffer

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In addition to the correct use and maintenance of the electrode, the electrode buffer that is matched with it also needs to be used and stored correctly

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After the pH buffer is prepared, it should be packed in a glass bottle or polyethylene plastic bottle

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For alkaline buffers, such as pH 9.18, pH 10.01, pH 12.46, etc., should be stored in plastic bottles

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In the refrigerator low temperature (5 ~ 10 °C) can generally be used for about 1 month

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If you find that the buffer is cloudy, moldy, or precipitated, it can not be used anymore

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Before calibrating with buffer, leave at ambient temperature for 1 to 2 hours for use after temperature equilibration

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Buffer should not be poured back after use to avoid contaminating other unused buffers

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Changes in temperature or chemical reactions in the sample itself can cause changes in pH. Therefore, you should try to keep the temperature consistent and avoid chemical reactions.

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Why is the buffer cannot be used within the validity period?

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The expiration date of the buffer refers to the shelf life in the unopened state of use. Once opened and used, the buffer is more susceptible to deterioration due to the action of various molds in the air.

## PA\_2022001\_00

---

Note: The buffer that has been used should never be poured back into the original bottle!

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How often do I need to calibrate my electrodes?

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The frequency of calibration of the electrode depends on the use, maintenance, nature of the sample, and measurement accuracy of the electrode. It is recommended to calibrate once a day, and do not exceed one calibration per week at the longest. Replacement of electrodes and prolonged non-use must be calibrated before use.

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The sample temperature is 10°C, at this time the meter displays the pH value at 10°C or 25°C?

The acidity meter shows the pH of the solution at the current temperature, if measured at 10 °C, the instrument shows the value of the solution 10 °C, if you need to get a pH of 25 °C, you must raise/cool the solution temperature to 25 °C, and then measure. Temperature compensation for acidity meters refers to compensating for the effect of temperature on the pH electrode, but cannot compensate the pH value at any temperature to 25°C.

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Sartorius

Copyright notice:

Sartorius reserves the right to make changes to the Technical Note without notice.

This Technical Note, including all of its components, is protected by copyright. Any use beyond the limits of the copyright law is not permitted without our approval. This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Sartorius