操作说明 | Operating Instructions

PB-30 / PB-30L

pH（酸度）计 | pH Meter

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This table was developed according to the provisions of SJ/T 11364.
本表依据 SJ/T 11364 的规定编制

X: Indicates that contents of Pb, Hg, Cd, Cr⁶⁺, PBB, PBDE in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T26572.
X: 表示该有毒或有害物质至少在该部件所用的某一均质材料中的含量超出 GB/T26572 标准规定的限量要求。

0: Indicates that contents of Pb, Hg, Cd, Cr⁶⁺, PBB, PBDE in at least one of the homogeneous materials used for this part is below the limit requirement in GB/T26572
O: 表示该有毒或有害物质在该部件所有均质材料中的含量均在 GB/T26572 标准规定的限量要求以下。

Content

1 About these Instructions	4	12.2 Package pH meter	27
1.1 Scope.....	4	13 Maintenance and protection	28
1.2 Symbols Used	4	13.1 Service	28
2 Safety Instructions	5	13.2 Clean pH Meter.....	28
2.1 Guidelines and basic information	5	14 Discard	29
2.2 Install information.....	5	15 Technical specifications	30
3 Device Description	7	15.1 General data	30
4 Getting Started	8	15.2 Specific model parameters.....	30
4.1 Unpacking and Equipment Supplied	8	15.3 Optional list of accessories.....	31
4.2 Mounting pH meter.....	8	15.4 Dimensions	32
4.3 Power connection	9	16 Declaration of Conformity	33
5 Use of pH meter	11		
5.1 Main interface operation display elements.....	11		
5.2 Temperature display and input.....	11		
5.3 Units of measurement settings	11		
6 Standardization and calibration	12		
6.1 Input electrode number	12		
6.2 Select calibration buffer	12		
6.3 Calibration	12		
6.4 Calibration report interface	12		
7 Data management	14		
7.1 Historical calibration data.....	14		
7.2 Historical measurement data	14		
7.2.1 Print historical data range.....	15		
8 User settings	16		
8.1 Set the user for the first time.....	16		
8.2 Create	16		
8.3 Edit	17		
8.4 User administration	18		
8.5 PB-30L password enable function	19		
9 System setting	20		
10 Communication interface	21		
10.1 USB-A port.....	21		
10.1.1 Scanner	21		
10.1.2 U-disk.....	21		
10.2 USB-B port.....	22		
10.2.1 Date output to Excel	22		
10.2.2 PC-SBI	23		
10.3 DB25 port.....	23		
10.4 SBI communication protocol	24		
10.4.1 Data input.....	24		
11 ISO/GLP compliant printout	26		
12 Transport pH Meter	27		
12.1 Remove pH Meter	27		

1 About these Instructions

1.1 Scope

These instructions are part of the device. These instructions apply to the device in the following versions:

Device	Model name
pH Meter, with or without User administration, electrode serialNo., USB-A port Function	PB-30 PB-30L

1.2 Symbols Used

DANGER

Denotes an immediate danger with risk that death or severe injury will result if it is **not** avoided.

WARNING

Denotes a danger with risk that death or severe injury may result if it is **not** avoided.

NOTICE

Denotes a danger that can result in property damage if the risk is **not** avoided.

Figures on the Operating Display

The figures on the operating display of the device may differ from those described in these instructions.

2 Safety Instructions

2.1 Guidelines and basic information

- This pH meter complies with the relevant directives and standards regarding electrical safety and electromagnetic compatibility. However, improper use or handling of equipment may result in damage and/or damage. Any improper use or operation of the pH meter, i.e. failure to use or operate in accordance with the manual, will result in the loss of all claims under the manufacturer's warranty.
 - Users need to read and understand these installation instructions, including safety instructions.
 - When used in systems and environments with higher safety requirements, you must comply with the requirements and regulations of your country/region.
 - Always keep your device and pH meter in a place where they are convenient to use.
- Any installation work or pH meter operation that fails to comply with the instructions will result in the loss of all claims under the manufacturer's warranty.

Risk of explosion!

Do not use the equipment in dangerous areas with explosive materials.

NOTICE

Do not use the equipment in dangerous areas with explosive materials.

2.2 Install information

WARNING

If the case or AC adaptor is damaged, do not operate the pH meter. Pull the plug out immediately to disconnect the damaged device from the power supply.

NOTICE

Do not expose the pH meter, AC adaptor of the pH meter or accessories coming with Sartorius to extreme temperature, corrosive chemical vapor, moist, impact, shock, or strong electromagnetic fields. Please observe the operating conditions described in the technical specifications.

NOTICE

The operator shall be responsible for any modification of the equipment and any connection of cables or equipment not provided by the Sartorius. Sartorius will ask for information about operational quality. You should only use peripherals provided by Sartorius.

Please note the IP protection level of the pH meter and its AC adaptor. Please do not allow liquids to penetrate. The level of protection indicates the adaptability of the equipment under various environmental conditions (humidity, foreign matter).

WARNING

If the case or AC adaptor is damaged, do not operate the pH meter. Pull the plug out immediately to disconnect the damaged device from the power supply.

The pH meter can only be opened by trained professionals of Sartorius.


Do not open the AC adaptor.

Risks during installation and operation

NOTICE

When laying the cable, be careful to choose the right position to prevent people tripping over the power line.

3 Device Description

Pos.	Name and Function	Pos.	Name and Function
1	Electrode support	6	Power connector
2	Electrode support base	7	Electrode support base
3	Touch screen	8	Temperature compensation joint
4	USB A port*	9	DB-25
5	Electrode	10	USB B port

Note: * means PB-30L does not contain this function.

4 Getting Started

4.1 Unpacking and Equipment Supplied

Open the packaging, making sure to remove all parts carefully.

After unpacking the device, check it immediately for any external damage.

If you detect any damage, proceed as directed in the chapter entitled "Care and maintenance."


Save the box and all parts of the packaging for any future transport. During shipment, please do not leave cables plugged in!

The following parts are included in the equipment supplied:


Components	PB-30 PB-30L
Host	1
Electrode support base	1
Electrode support	1
Electrode	1
* Complete set of pH buffering agent	2 set
AC adaptor with power connector specific to each country	1
User manual	1

Note: * Optional according to transportation regulations of different countries.

4.2 Mounting pH meter

Mount the support on the base.

Mount the electrode on the support and the pH meter mounting is finished.

4.3 Power connection

WARNING

Use of an incorrect power plug adaptor or improper use of the power plug adaptor may result in a fatal electrical shock.


Connect the power plug adaptor to the AC adaptor specific to a country/region. The power plug adaptor must be suitable for wall-mounted power sockets at the installation location.

Do not plug the power adaptor into the socket without the AC adaptor.

Connect the pH meter to the AC power supply

Use only the original Sartorius AC adaptor.

For use within China only: Part No. 6971793S


- 1) First, plug the AC adaptor into the power socket on the back of the pH meter.
- 2) Next, connect the AC adaptor to the wall-mounted socket (main power supply).

Power connection (backup power)

Item No. on the package	Power supply/country-specific power plug adaptor (packaged in the polyethylene bag printed with country codes such as EU)	Illustration (from left to right)
YEPS01-PS4	Power supply with connecting cables	
	China (CN) South Africa (ZA) Australia (AU)	
	India (IN) The United States and Japan (US+JP) Europe (EU)	
	The United Kingdom (UK)	
YEPS01-PS5	Power supply with connecting cables	
	Argentina (AR) Brazil (BR) Korea (KR)*	

Note: * means that YEPS01-PS4 European (EU) connector can be used instead of Korean products.


Select power plug adaptors for specific countries / regions. The power plug adaptor must be suitable for wall-mounted power sockets at the installation location.

Push the power plug adaptor into the fixator of the AC adaptor. The ribbed key must be facing forward.

Push the power plug adaptor in completely until you hear the insertion sound.

Check if the power plug adaptor is firmly locked in place. To this end, please gently pull the power plug adaptor.


If the power plug adaptor cannot be moved, it indicates that the adaptor has been locked in place.

Remove the power plug adaptor

Press the ribbed key from above and slide back into the power plug adaptor.

Push out and remove the power plug adaptor from the AC adaptor.

Connect the AC adaptor

Check the rated voltage on the nameplate of the AC adaptor. Ensure that the voltage level printed on the device meets the local supply voltage at the installation location.

If the calibrated supply voltage fails to meet the local supply voltage or there is no suitable AC adaptor available:

Do not use the AC adaptor. Please contact Sartorius Service.

Use only the original Sartorius AC adaptor.

Connect the angle plug (2) to the power socket of the device (operating voltage connection).

Connect the power cord (1) to the wall-mounted socket (power voltage) in the mounting position.


5 Use of pH meter

This section describes the options you can use to familiarize yourself with the pH meter.

5.1 Main interface operation display elements

NOTICE

Sharp objects (such as ballpoint pens) may cause damage to the equipment!

1	User time bar*	8	Standard buffer (relative mv offset)
2	Temperature mode	9	Calibration
3	Electrode serial No.*	10	GLP header and footer
4	Timing calibration	11	Print
5	Temperature display	12	Measured value
6	Data management	13	Units of Measurement (pH, mV, RmV)
7	System setting	14	Alarm

Note: *PB-30L does not display user name or electrode serial number.

5.2 Temperature display and input

The screenshot shows the 'Temp.compensation' screen. At the top, there is a status bar with 'User:', '2017-05-03', and '12:34:58'. Below this, there are buttons for 'Cancel' and 'Enter'. The main display area shows 'TempCM Value:' followed by an input field containing '60.0' and a range indicator '(-5.0°C~105.0°C)'.

When the device is connected to the temperature sensor, the device enters the automatic temperature compensation state, and the temperature compensation mode is displayed as ATC;

When the device is not connected to the temperature sensor, the device enters the manual temperature compensation state, and the temperature compensation mode is displayed as MTC. At this point, click the temperature indication display area of the main interface, you can enter the manual temperature compensation page, enter the temperature compensation value, and confirm.

5.3 Units of measurement settings

The device can display three measurement units, pH, mV, RmV (relative mV offset).

You can switch the measurement unit by clicking the measurement unit display area on the main interface.

Any measurement unit can be turned on/off through the unit settings in the system setup menu.

6 Standardization and calibration

Before using the pH meter, it shall be calibrated to ensure accurate measurement data.

6.1 Input electrode number

When the Units of Measurement is pH, click the "calibration" key on the home screen and the system enters the electrode No. input interface*.

Note: * means PB-30L carries no such function.

6.2 Select calibration buffer

Select a set of calibration buffers, which should agree with the pH value of the calibration solution, otherwise there will be a calibration error, or there will be measurement data errors after calibration.

6.3 Calibration

The system carries out calibration of up to three points. When to calibrate, put the electrode in the calibration solution and the pH meter will automatically identify the current value of the calibration buffer, which is blinking black and white; click "Confirm" to start the calibration, until the value no longer flashes, indicating that the calibration is finished.

After calibration of one point at least, a key **Finish** will appear on the interface. Click the key and the system will enter the calibration report interface.

Click on the unit of measurement to switch between "pH" and "mV". If the unit of measurement is "mV", the key "Confirm" will turn grey. Only when the unit is pH can the key "Confirm" be clicked for calibration.

Click the key "Cancel" to have the system enter the calibration report interface.

6.4 Calibration report interface

User:

2017-05-03

12:34:58

Calibration report

Calibration report 3:

Std. buffer

:

S2

25.0°C

ATC

6.86pH

5.0mV

Std. buffer

:

S3

25.0°C

ATC

9.18pH

-305.0mV

←

^

v

📄

User:

2017-05-03

12:34:58

Calibration report

Calibration report 3:

Electrode slope

:

Slope1

97.6%

Electrode slope

:

Slope2

95.3%

Calibration status

:

Calibration OK


←

^

📄

Click  to save the calibration data and return to the home screen.



7 Data management

The pH meter model PB-30 can record 30 sets of historical calibration data, 400 sets of historical measurement data, and 50 sets of log data.


The acidity meter model PB-30L can record 10 sets of historical calibration data, 500 sets of historical measurement data, no log data.

7.1 Historical calibration data


Click and to view different sets of historical reports;

click and to view the current set record information.


7.2 Historical measurement data

Click the "Historical measurement data" display bar to enter the historical measurement data interface.

Page up 50 sets of data.

Page down 50 sets of data.

Click on each piece of data and enter the "Historical data details" page.

7.2.1 Print historical data range

User: 2017-05-03 12:34:58

< Data Measurement report

Print item range: (eg: 1-10)

1	2	3	←	✕
4	5	6		
7	8	9	C	✓
0	-			

Press the print key to display this page.

Print item range: (eg: 1-10)

The data display box shows the range of printed data input by the user, requiring that the data following "-" must be larger than the previous data.

[-] This key is used to indicate the range of data from ... to ..., for example, 10-100 indicates printing historical measurement data from Set 10 to 100.

✕ Cancel print.

✓ Start print.

7.3 Log report

User: Admin 2017-05-03 12:34:58

< Date Log data

Log report 3:

Date and time : 2017-05-01 08:00:52
 User : Admin
 Content : Change language

< >

📄

Click < and > to scroll through the historical report information for different groups.

Click 📄 to save all log reports in a PDF file on the USB flash drive at one time, and click again to save all log reports in a new PDF file.

8 User settings

PB-30 user settings allow setting user name, user permissions and password functions, etc., PB-30L user settings only allows setting whether password is enabled.

A maximum of 5 users can be set up, and there must be an administrator and a tourist. The administrator name is fixed as Admin and the tourist name is fixed as Guest. You cannot modify the name. When creating an Admin, a guest is automatically created.

Guest users can only perform measurements and cannot perform other operations.

Enter the user setting mode through the main page Setting.

8.1 Set the user for the first time

When setting a user for the first time, the user list is empty

Click **Create** to create a new user. Since it is the first time to create a new user, the user must be an administrator.

Click **Edit** to edit the user information, such as modifying the user name, changing the password, and deleting the password.

Modify cannot be used when there is no user.

When there is no user, click to create, and only the administrator can be created. The administrator's name can only be Admin.

Click "Next" to enter the password enable interface.

With password disabled, click "Next" to complete the user creation.

Enable the password and click "Next" to enter the password setting interface.

The user password shall be less than 14 digits, which can be a combination of numbers and letters.

8.2 Create

Only user Admin can create new users. Click **Create** to create a new user.

Click "Next" to confirm the use of the current user name and enter the user password enter interface;

Click the "USER-1" display bar to enter the "Enter Username" interface and enter the user name.

User name can only be composed of letters + numbers, up to 14 digits.

Set user right.

Click "Create 2" to return to the "Enter Username" interface.

With password disabled, click "Next" to complete the user creation.

Enable the password and click "Next" to enter the password setting interface.

The user password shall be less than 14 digits, which can be a combination of numbers and letters.

8.3 Edit

Before the user enters the user edit interface, if there is a password, the user shall enter the password first correctly before entering the user edit interface.

User edit allows modifying the user name, setting password and deleting password.

* Admin cannot modify its own username.

User: User1	2017-05-03	12:34:58
< User list Edit		
Edit user name >		
Edit user password >		
Delete user password >		

8.4 User administration

Only when the user is an administrator, can Admin manage the users. Before entering user administration and if Admin has a password, it is required to enter the password before entering the user administration interface.

User: Admin	2017-05-03	12:34:58
< User list Modify		
Modify user rights >		
Delete user >		
Delete user password >		

User: Admin	2017-05-03	12:34:58
< Cancel	User rights Next	
User 1	✓	
User 2		
User 3		

Select a user and click Next to modify user permissions.

User: Admin	2017-05-03	12:34:58
< Modify	User rights Enter	
Language	🔘	
Date and Time	⬆	🔘
Units	▼	🔘
Timing Cal.	🔘	

Modify user rights.

User: Admin	2017-05-03	12:34:58
< Modify	User rights Enter	
Printout	🔘	
Interfaces	⬆	🔘
Calibration	⬇	🔘

User: Admin	2017-05-03	12:34:58
Cancel	Delete user	Enter
Admin		✓
User1		✓

If the administrator is deleted, all users will be deleted by default.

User: Admin	2017-05-03	12:34:58
Delete User		
Do you want to delete user? Admin, User1		
✓ ✕		

User: Admin	2017-05-03	12:34:58
Cancel	Delete user password	Enter
Admin		✓
User1		✓

Select each user, and the user will be followed by a checkmark, indicating that the user is selected. Click the selected user again, and the checkmark will disappear, indicating that the user is unselected.

User: Admin	2017-05-03	12:34:58
Delete user password		
Do you want to delete user password ? Admin, User1		
✓ ✕		

8.5 PB-30L password enable function

PB-30L supports no other functions of user settings, other than password enable.


User:		
< Setting	Password	Next
Yes, enable password		✓
No, disable password		

To enable password protection, you need to enter the password when you click System Settings.

User:		
< Setting	Password	Enter
Yes, enable password		
No, disable password		✓


Password protection is not enabled. You do not need to enter a password to click system settings.

9 System setting

Note: *PB-30L carries no such function.

10 Communication interface

When the equipment is connected to a computer, the operating systems that are required include Windows 7, Vista, Windows XP, or 2000.

User:	2017-05-03	12:34:58
< Setting Interfaces		
PC (USB B)	output to Excel	☒
Printer (COM)	YDP20	☐
Peripheral interface (USB A)	Scanner	☐

The default is that the USB B port is enabled to have the data exported to Excel, and the other ports are disabled by default.

Note: PB-30L carries no peripheral interface (USB A port).

User:	2017-05-03	12:34:58
Cancel Peripheral... Enter		
Scanner		☒
U-disk		☐

10.1 USB-A port

USB-A can only be connected to one device at a time. The device can be used only if an appropriate function is selected. If the device is connected, but the port function is not enabled, the device cannot function properly.

10.1.1 Scanner

The scanner can identify barcode and QR code, and it can be directly used by being connected to the device with no need of driver.

The scanner can replace the keyboard for entering electrode serial number and sample ID.


10.1.2 U-disk

The name of the file in the USB-disk is made up of the date and time of creation, for example: "2018-06-06 12-13-36.txt".

Through the print function, the data will be automatically stored on the USB-disk in the printing format as required.

The USB flash drive used is less than 8GB in storage capacity.

Example:

The figure contains GLP-printing header and footer, as well as historical calibration data, and measurement data.

10.2 USB-B port

The USB-B port connects to the PC and can output data to an Excel spreadsheet. It can also communicate using the SBI protocol.

10.2.1 Date output to Excel

With this feature, the device will be virtual to the computer's keyboard peripheral to transfer data to the computer. After the device is connected to the computer, open an Excel table, and it is recommended to set the font to Courier New to facilitate data alignment; move the cursor to the location you want to output, switch the input method into English, and click on the print key on the screen, the data will be automatically uploaded to the location where the cursor is in the table.

User:	2017-05-03	12:34:58
Candle	PC	Enter
Output to Excel		✓
USB to RS232 SBI		

Example:

	A	B	C
1	-----		
2	2018-05-17	15:45:55	
3	Sartorius		
4	Model	PB-30	
5	SerNo.	1234SERIAL1234	
6	VerNo.	00.01	
7	-----		
8	EtdNo.	123ELECTRODE34	
9			
10	2018-05-09	21:44:44	
11	User	User-000000004	
12	EtdNo.	SAR-0000-00004	
13	Calibration Clear		
14	2018-05-08	09:43:00	
15	User	User-000000003	
16	EtdNo.	SAR-0000-00003	
17	S1	25.2 °C	
18	ATC		
19	7.02 pH		
20	0.5 mV		
21	Slope1	100.0 %	
22	Calibration OK		
23			
24	-----		
25	2018-05-17	15:46:31	
26	Name:		
27			
28	-----		

The table shows the history calibration information with GLP header and footer.

10.2.2 PC-SBI

With this function, the device will act as the RS232 interface and automatically match the RS232 interface settings to communicate with other devices. After the device is connected to the computer, the computer will automatically search for the driver and install it, so that the user can use it directly. If the driver installation fails, please go to the website ... to download the driver and install it on the computer.

Open the serial port for data transmission.

The data transmission content complies with the printing requirements and displays the information that the pH meter needs to print.

10.3 DB25 port

DB25 port allows connecting to two different types of printers for data print. And it can be used as an RS232 interface device to communicate with other RS32 interface devices through the SBI communication protocol.

The port is set to by default: Baud rate: 9600

Data bits: 7

Stop bits: 1

Check bits: odd

User:	2017-05-03	12:34:58
Cancel	Printer	Enter
YDP20...		✓
YDP30		
RS232 SBI		

The terminal DB25 is connected to the pH meter; the terminal DB9 is connected to the printer or the computer.

10.4 SBI communication protocol

Users can read the equipment information and real-time measurement data of the pH meter by upper computer software according to the SBI communication protocol.

10.4.1 Data input

Interface command

The computer connected through the data port sends control commands to the pH meter to control its function.

Format 1	Esc	!	CR	LF		
Format 2	Esc	!	#	_	CR	LF

Esc	Exit
!	Command character
#	Figure
_	Underline (ASCII:95)
CR	Enter (optional)
LF	Line feed (optional)

Example:

Real-time measurement data read command:

ESC P CR LF

Regular data return format

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
K	K	K	K	K	K	+	A	A	A	A	A	A	A	A	A	*	E	E	E	CR	LF
K	K	K	K	K	K	*	A	A	A	A	A	A	A	A	A	A	A	A	A	CR	LF
A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	A	CR	LF

K	Mark
CR	Carriage return
*	Space
LF	Line feed

A Display character
E Measurement unit mark

Example:

1. Currently measured value by user:

Value 10.02 pH

2. Instrument serial number:

SerNo XXXXXXXXX

Command format (character form)	Function
ESC x0_ CR LF	Read the instrument manufacturer
ESC x1_ CR LF	Read the instrument model
ESC x2_ CR LF	Read the instrument serial number
ESC x3_ CR LF	Read the instrument software version number
ESC P CR LF	Read real time measurement data
ESC P1 CR LF	Read date and time data
ESC P2 CR LF	Read electrode ID
ESC P3 CR LF	Read sample ID
ESC P4 CR LF	Read the temperature compensation mode
ESC P5 CR LF	Read the temperature
ESC P6 CR LF	Read the pH value
ESC P7 CR LF	Read the mV value
ESC P8 CR LF	Read the RmV value
ESC P9 CR LF	Read the user ID
ESC Q CR LF	Read all historical measurement data
ESC R CR LF	Read all historical calibration data

11 ISO/GLP compliant printout

The device data, ID and current information can be printed before (GLP header) and after (GLP footer) the value of the measurement series component.

GLP header printing contains:

1. The current printed date and time
2. Sartorius logo
3. Product model
4. Product SerNo.
5. Software VerNo.
6. User name
7. Electrode No.

GLP footer printing contains:

Blank line " "

1. The current printed date and time
2. Signature block "Name:"

Print example:

```
-----
2018-05-17  15:45:55
      Sartorius
Model          PB-30
SerNo.1234SERIAL1234
VerNo.          00.01
-----
EtdNo.123ELECTRODE34
```

```
-----
2018-05-17  15:46:31
Name:
```

```
-----
```

```
-----
2018-05-17  15:45:55
      Sartorius
Model          PB-30
SerNo.1234SERIAL1234
VerNo.          00.01
-----
EtdNo.123ELECTRODE34
```


```
2018-05-06  21:41:00
User  User-000000001
EtdNo.SAR-0000-00001
S1          25.0 °C
           ATC
           7.00 pH
           0.3 mV
Slope1      100 %
Calibration OK
```

```
-----
2018-05-17  15:46:31
Name:
```


```
-----
```

12 Transport pH Meter

12.1 Remove pH Meter

First remove the electrode from the support


Next, remove the support from the base in the end.


12.2 Package pH meter

Step 1: put the pH meter into the foam package;

Step 2: fold the left and right sides up and wrap the host completely;

Step 3: put the packaged motor core into the packing box, then put the support, base, electrode and power supply into the packaging box, and close the box cover.

13 Maintenance and protection

13.1 Service

To ensure a stable measurement accuracy of the pH Meter, we recommend that you arrange regular maintenance at least once a year. The AC adaptor and its connection must be checked regularly by a professional electrician (for example, every 2 years).

Maintenance work can only be performed by trained maintenance technicians. The plug of equipment must be pulled out during maintenance and it is required to pull the plug out of the main power supply. Besides, the accuracy of the pH Meter is likely to be affected and may pose a serious risk to users. Please contact the Sartorius service or the Sartorius dealer for proper maintenance of the pH Meter.

13.2 Clean pH Meter

WARNING

Contaminated equipment:

- Health risks posed by the existence and collection of residual microbial contamination in the product.
- Health risks posed by biological or microbial substances.
- Please observe the cleaning technical specifications.
- Please check the cleaning results carefully.

WARNING

- Disconnect the AC power from the device: unplug the AC adaptor from the power supply.
- Disconnect all data cables connected to the pH meter.

NOTICE

- Do not open the pH meter or AC adaptor.


The parts contained in these devices cannot be cleaned, maintained or replaced by the operator.

- Make sure no liquid or dust enters the pH meter or AC adaptor.
- Do not use detergents containing solvents or abrasive ingredients, which may eventually cause damage to the equipment.

NOTICE

Do not use acetone or corrosive cleaning agents to clean the following parts:
Control panel, power plug, data interface, label, and all other plastic parts.

14 Discard

The outer packaging is made from environmentally friendly materials and be used as secondary raw materials. In Germany, the packaging out of use can be handled for free through the VFW dual system (contract No. D-59101-2009-1129). In other countries or regions, please use local waste disposal facilities to dispose of materials in accordance with applicable regulations. Equipment including accessories and batteries shall not be treated as general household waste and they shall be recycled as electrical and electronic equipment. For disposal or recycling, please contact our field service personnel. Besides, in addition, the partners listed in the following websites can deal with the waste within the EU:

- 1) Select <http://www.sartorius.com>.
- 2) Select the "Service" tab.
- 3) Then select "Disposal information".
- 4) Download the PDF file on this page, which contains the information of local Sartorius addresses for disposal.

Sartorius does not recover equipment contaminated by hazardous substances (ABC pollution), nor provides maintenance or disposal services.

Disposal service address

For more information about the service address or equipment disposal, please visit our website (www.sartorius.com) or contact the Sartorius service department.

15 Technical specifications

15.1 General data

Specification	Unit	Value
AC adaptor		
Sartorius AC adaptor		
AC input		100-240 V ~ , -10%/+10% , 50-60 Hz , 0.2 A.
DC output		15 V DC , ±5%, 530 mA (max) / 8 W (max): 0 to +40 °C ; 15 V DC , ±5%, 330 mA (max) / 5 W (max): 0 to +50 °C .
Other Data		Protection level II EN/IEC 60950-1 / over sea level 3000m at most ; IP40 EN 60529/IEC 60529.
pH meter		
Input voltage	V DC	12.0-15.0
Power consumption	W	4.5 (general)
Environmental conditions		
This specification is used under the following environmental conditions:		
Environment		For indoor use only.
Ambient temperature	°C	+10 - +30
Run range	°C	Guaranteed between +5 and +45.
Storage and shipping	°C	-10 - +60
Height	Altitude (m)	Up to 3000
Relative humidity		When the temperature is below 31°C, the relative humidity is 15-80%, without condensation; when the temperature reaches 40°C, the relative humidity decreases linearly to 50%; when the temperature reaches 50°C, the relative humidity is 20%.
Safety of electrical apparatus		In accordance with EN 61010-1/IEC 61010-1 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use Part 1: General Requirements.
Electromagnetic Compatibility		In accordance with EN 61326-1 / IEC61326-1 Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements - Part 1: General Requirements.
Anti-interference		Applicable to industrial areas.
Transient emission		Level B, suitable for residential areas or areas of low voltage network for residential buildings.
Standard equipment		
Measurement units optional		pH, mV and RmV


15.2 Specific model parameters

Measurement parameters	pH	mV
Measurement capability	-1.99 - 19.99 pH	-1999.9 - 1999.9 mV
Measurement accuracy	± 0.005 pH	± 0.3 mV
Read accuracy	0.01pH	0.1mV

15.3 Optional list of accessories

Electrode	As need
Calibration solution	As need
Support	As need
Base	As need
Printer	YDP20-OCEV1 Dot matrix printer
Printer	YDP30 Thermal/thermal transfer printer
YDP30 printing paper	69Y03285
USA - A port concerted to USB - B port	As need
DB25 converted to DB9 cable	YCC03 - D25
Power adaptor	6971793S / YEPS01 - PS4 / YEPS01 - PS5
Scanner	As needed

15.4 Dimensions

All dimensions are given in millimeters.

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in the form as delivered fulfils all the relevant provisions of the following European Directives – including any amendments valid at the time this declaration was signed – and meets the applicable requirements of the harmonized European Standards including any amendments valid at the time this declaration was signed – listed below:

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The person authorised to compile the technical file:

Sartorius Lab Instruments GmbH & Co. KG
Electronics & Product Compliance
37070 Goettingen, Germany

Sartorius Scientific Instruments (Beijing) Co., Limited
Beijing, 2019-05-23

Dr. Reinhard Baumfalk
Head of Product Development
Lab Products and Services Division

Zhu Xiulin
R&D Engineer, International Certification Management

Sartorius Lab Instruments GmbH & Co. KG
Otto-Brenner-Strasse 20
37079 Goettingen, Germany
Phone +49.551.308.0
Fax +49.551.308.3289
www.sartorius.com

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The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Copyright notice:

This instruction manual, including all of its components, is protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:
12|2019

Printed in the CN on paper bleached
without chlorine. | JB
Publication No.: WPB6005-aw191202